

Select the instruments you need to perform a complete diagnostic examination:



Green Series™ 777 Wall Transformer
Energy efficient with handle-based rheostats and automatic on/off sensors make operation simple.



3.5 V Coaxial Ophthalmoscope
Multiple aperture/filter options provide long-lasting, reliable performance for general and specialist examinations.




PanOptic™ Ophthalmoscope
5x larger field of view of the fundus vs. standard ophthalmoscopes in an undilated eye (25° vs. 5°).



MacroView™ Otoscope
A nearly complete view of the tympanic membrane with 30% greater magnification than a traditional otoscope.



Pneumatic Otoscope
Rotating lens with ultraseal enables pneumatic otoscopy and easier instrumentation.



KleenSpec® Disposable Specula Dispenser
Offers optimum organization of disposable otoscope tips. Comes filled with specula.



Wall Aneroid
Easy-to-read, jeweled-movement, mercury-free manometer is certified to be accurate within one percent of scale (± 3 mmHg).



**Spot Vital Signs®/Spot Vital Signs® LXi/
Connex® ProBP™ Devices**
Fast, accurate, and efficient way to take vitals. Allows for electronic transfer of data to EHR.




SureTemp® Plus Electronic Thermometer
Provides accurate oral readings in about four to six seconds, 10- to 13-second pediatric axillary and rectal readings and 12- to 15-second adult axillary readings.

The Welch Allyn Green Series™ 777 Integrated Wall System

WelchAllyn®

Advancing Frontline Care™

Recommended Systems:



77791-MX
Includes GS 777 Wall Transformer (#77710), Coaxial Ophthalmoscope (#11720), MacroView Otoscope (#23810), Wall Aneroid (#7670-01), KleenSpec Dispenser (#52400-PF), and SureTemp Plus Thermometer (#01690-300).



77791-2MPX
Includes GS 777 Wall Transformer (#77710), PanOptic Ophthalmoscope (#11820), MacroView Otoscope (#23810), Wall Aneroid (#7670-01), KleenSpec Dispenser (#52400-PF), and SureTemp Plus Thermometer (#01690-300).



77791-2MPXPROBP
Includes GS 777 Wall Transformer (#77710), PanOptic Ophthalmoscope (#11820), MacroView Otoscope (#23810), KleenSpec Dispenser (#52400-PF), Connex ProBP Digital Blood Pressure Device (#34XFWT-B) and SureTemp Plus Thermometer (#01690-300).



77794-2MPLXI
Includes GS 777 Wall Transformer (#77710), PanOptic Ophthalmoscope (#11820), MacroView Otoscope (#23810), KleenSpec Dispenser (#52400-PF), and Wall Mount for Spot Vital Signs LXi (#4500-62). Spot Vital Signs/Spot Vital Signs LXi ordered separately.

	WALL BOARD	GS 777 WALL TRANSFORMER	COAXIAL OPHTHALMOSCOPES	PANOPTIC™ OPHTHALMOSCOPES	DIAGNOSTIC MACROVIEW™ OTOSCOPES	PNEUMATIC OTOSCOPE	KLEENSPEC® DISPENSERS	WALL ANEROID	CONNEX® PROBP™ DEVICE	SURETEMP® PLUS THERMOMETER
Part #		77710	11720	11810 11820	23810 23820	20200	52400-PF 52100-PF	7670-01	34XFWT-B	01690-300
GS 777 Integrated Wall Systems										
77798	•	•	•			•		•	•	
77791-1MPX	•	•		•	•		•	•		•
77791-2MPX	•	•		•	•		•	•		•
77791-MX	•	•	•		•		•	•		•
77792-1MP	•	•		•	•		•	•		
77792-2MP	•	•		•	•		•	•		
77792-M	•	•	•		•		•	•		
77795-2MPX	•	•		•	•		•	•		•
77795-MX	•	•	•		•		•	•		•
77796-M	•	•	•		•		•	•		
77797-2PX	•	•		•		•		•		•
77797-X	•	•	•			•		•		•
77791-MXNOBP	•	•	•		•		•			•
77792-MNOBP	•	•	•		•		•			
GS 777 Integrated Wall Systems with Connex® ProBP™ Digital Blood Pressure Device (includes wall mounting bracket and device)										
77791-MXPROBP	•	•	•		•		•		•	•
77791-2MPXPROBP	•	•		•	•		•		•	•
GS 777 Integrated Wall Systems for Spot Vital Signs® or Spot Vital Signs® LXi Devices (includes device wall mounting bracket only, device(s) must be ordered separately)										
77794-MSPOT	•	•	•		•		•			
77794-MLXI	•	•	•		•		•			
77794-2MPSPOT	•	•		•	•		•			
77794-2MPLXI	•	•		•	•		•			
GS 777 Wall Transformer Sets										
77510		•					•			
77910		•					•	•		
77710-71M		•	•		•					
77710-81M		•		•	•					
77710-82M		•		•	•					
77710-TAKX		•					•	•		•