



360 Route 101, Bldg 3, Bedford, NH 03110  
Tel: (603) 471-1499 [info@omnigloves.com](mailto:info@omnigloves.com)

---

---

## OmniShield™ Product Information Bulletin

---

---



### Anti-Microbial Vinyl Powder Free GP Gloves

#### PRODUCT DESCRIPTION

ANTI-MICROBIAL

DOP (di-octyl phthalate) / DEHP (di-2-ethylhexyl phthalate) FREE

NON-LATEX, AMBIDEXTROUS

POWDER FREE VINYL GENERAL PURPOSE GLOVE

#### PRODUCT REORDER NUMBERS

Size	Reorder Number
Small	352-11
Medium	352-12
Large	352-13
X-Large	352-14

#### Technical Information

<b>MATERIAL</b>	Polyvinyl Chloride and compounding materials safe for use general purpose gloves, DOP (di-octyl phthalate) / DEHP (di-2-ethylhexyl phthalate) free.
<b>GLOVE DESIGN</b>	Ambidextrous, with beaded cuff. 9.5 inch length.
<b>TESTING METHODS</b>	ASTM D5250-00
<b>STANDARDS</b>	ISO 2859
<b>ANTIMICROBIAL</b>	<b>AMz4.0</b>

Technical Information (continued)  
 Anti-Microbial Vinyl Powder Free GP Glove

<b>Packaging</b>	100 gloves per box, 10 boxes per case.  Dispenser Box Dimensions: 235mm x 127mm x 60mm Approximate weight: 1.5 lbs  Outer Carton Dimensions: 315mm x 260mm x 246mm Approximate weight: 15 lbs	
<b>Minimum Thickness(Single Thickness)</b>	Finger Tip Palm Cuff	0.09mm ± 0.03 0.08mm ± 0.03 0.08mm ± 0.03
<b>Minimum Length</b>	Small Medium Large X-Large	240 mm ± 10 240 mm ± 10 240 mm ± 10 240 mm ± 10
<b>Minimum Width at Palm</b>	Small Medium Large X-Large	85mm ± 5 95mm ± 5 105mm ± 5 115mm ± 5
<b>Estimated Glove Weight</b>	Small Medium Large X-Large	4.0 ± 0.5 grams 4.5 ± 0.5 grams 5.0 ± 0.5 grams 5.5 ± 0.5 grams
<b>Tensile Strength</b>	<u>Unaged</u> ≥9 (MPa)	<u>Aged</u> 9 (MPa)
<b>Elongation</b>	<u>Unaged</u> ≥300%	<u>Aged</u> ≥300%
<b>Residual Protein Content ASTM D3578 &amp; D5712 Standard</b>	Contains No Natural Rubber Latex	
<b>Powder Content</b>	<2mg total powder per glove	
<b>AQL</b>	4.0	